

## APPLICATION AREAS: ELECTRODES

APPLICATION	GE 100	GE 101	GE 104	GE 108	GE 114	GE 117	GE 120	GE 125	GE 151	GE 171	GE 173	GR 105	GR 175
Sewage											•		
Aquarium water	•		•	•	•	•			•			•	•
Soil testing		•											
Emulsions		•	•										
On-site measurements				•	•	•		•				•	
Fish farming	•		•	•	•	•		•	•			•	•
Galvanic baths											•		•
Beverages								•	•		•	•	•
Low-ion media			•								•		
Cosmetics			•										
Food sample		•					•						
Sea water	•	•	•	•	•	•	•	•	•	•	•	•	•
Online measuring										•	•		•
Process chemistry									•	•	•		•
Swimming pool water	•			•	•	•		•			•	•	•
Suspensions		•	•										•
Drinking water	•		•	•	•	•		•			•	•	•
Water-based lacquers			•								•		

Note: The set information are to provide general recommendations. It needs to be checked, which electrodes for each area of application are suitable.