## **GMH3511 SETS**



COMFORTABLE MEASUREMENT

#### **GMH3531 SETS**



FOR LABORATORY AND FIELD

#### GMH3551 SETS



WITH DATA LOGGER

## **GMH 3511-SET**

Art. no. 605021

pH-/Redox-/Temperature measuring instrument Device complete with pH electrode GE 114, GF1T-T3-B-BS, 5x GPH4, 5x GPH7, 2x GPF100

#### GMH 3511-G125

Art. no. 475740

pH-/Redox-/Temperature measuring instrument Device complete with pH electrode GE 125 (PT1000)

#### General:

For comfortable measurement of pH value and temperature. Even easier operation with a menu reduced to 5 points in GMH 3511

Minimum measurement effort with maintenance-free gel electrodes and automatic temperature compensation.

# Specifications:

see GHM 3511

## Accessories and spare parts:

## **GE 114-BNC**

Art. no. 604701 pH-electrode

## GE 125-BNC

Art. no. 600731

waterproof pH-electrode incl Pt1000, BNC plug

## **GKK 3001**

Art. no. 611605

with cut-outs for 1 device of the GMH 3xxx series and accessories for water analysis (395 x 295 x 106 mm)

#### **GMH 3531-SET125**

Art. no. 474240

pH-/Redox-/Temperature measuring instrument Device complete GE 125, 5x GPH4, 5x GPH7, 2x GPF100, GKK 3001

## **GMH 3531-SET100**

Art. no. 604591

pH-/Redox-/Temperature measuring instrument Device complete GE 100, GF1T-T3-B-BS, 5x GPH4, 5x GPH7, 2x GPF100, GKK 3001

#### ---

Functional scope for elevated demands in the laboratory and field

Minimum measurement effort with maintenance-free gel electrodes and automatic temperature compensation.

# Specifications:

see GHM 3531

## Accessories and spare parts:

## GE 100-BNC

Art. no. 600704

pH-electrode, BNC plug

## GE 125-BNC

Art. no. 600731

waterproof pH-electrode incl Pt1000, BNC plug

## **GKK 3001**

Art. no. 611605

with cut-outs for 1 device of the GMH 3xxx series and accessories for water analysis (395 x 295 x 106 mm)

## **GMH 3551-SET125**

Art. no. 474903

pH-/Redox-/Temperature measuring instrument with logger; Device complete GE 125, 5x GPH4, 5x GPH7, 2x GPF100, GKK 3001

## **GMH 3551-SET100**

Art. no. 475742

pH-/Redox-/Temperature measuring instrument with logger; Device complete GE 100, GF1T-T3-B-BS, 5x GPH4, 5x GPH7, 2x GPF100, GKK 3001

#### General:

Sets with very good equipment with integrated data logger

#### Specifications:

see GHM 3551

#### Accessories and spare parts:

## GE 100-BNC

Art. no. 600704

pH-electrode, BNC plug

## **GE 125-BNC**

Art. no. 600731

waterproof pH-electrode incl Pt1000, BNC plug

## GKK 3001

Art. no. 611605

with cut-outs for 1 device of the GMH 3xxx series and accessories for water analysis (395 x 295 x 106 mm)