


LOGGER SOFTWARE



HIGHLIGHTS:

- Automatic reading
- Operation of the logger function
- Diagram display
- Export function

GSOFT 40K

Art. no. 601145

Windows-software for EASYLog and T-Logg with logger

General:

Software for starting, stopping and reading of EASYLog devices. The data can be visualised, saved and prepared for further processing.

Application:

Incoming goods control and cooling chain, test stand and laboratory equipment monitoring.

Specifications:

Program version:	Application with user interface
Data backup:	File (binary)
Export formats:	*.csv
Languages:	German English French Czech
Access control:	-
Remote access:	-
Alarms:	-
Devices:	EASYLog (via EASYBus level converter) T-Log (via USB 100)
Severl interfaces:	Can be used in succession
Recording interval:	Depending on the data logger
Live display:	No
Reading of data loggers:	Yes
System requirements:	Microsoft Windows 7 SP1 (32 or 64 Bit) Not executable with Windows RT, Windows 10 im S Modus, ARM or Intel Itanium based Windows systems
Scope of supply:	DVD, manual

SOFTWARE CONNECTION

	GMH 3000-DLL (windows library)	inray - inMOVE GPL-Router-Plug-in (www.inray.de) (not shown)	Protocol description (www.greisinger.de) (w/o picture)
EASYBus and EASYLog	•	•	•
TLogg	•		•
Start, stop, delete, read-out logger	•		
Program examples	Visual Studio, Excel VBA		
subject to costs	•	•	

EASYBUS.dll

Art. no. 609174

Windows function library for communication with EASYBus devices